



Scintillator tests for SpecMAT

J. A. Swartz

ACTAR TPC Collaboration Workshop
GANIL, Caen
19/11/2015

Contents

1. Introduction
2. Cubic scintillators with SiPMs
3. Cylindrical scintillators with PMTs + GET
4. SiPMs + GET
5. Summary + Outlook

Co-authors:

KU Leuven - IKS: R.Raabe, H. De Witte, F. Renzi, T. Marchi

University of York: D. Jenkins, N.S. Bondili, P. Joshi

GANIL: G.F. Grinyer

The SpecMAT project:

To investigate exotic nuclei, specifically

- neutron-rich around **Z = 28**, e.g. $^{68, 70}\text{Ni}$; $^{78, 80}\text{Zn}$ (d, p); ($d, ^3\text{He}$)
- neutron-deficient around **Z = 82**, e.g. ^{184}Hg , ^{188}Pb , ^{196}Po (d, p); (p, d)
 - To place active target detector inside solenoid magnet,
 - & supplement with scintillator array in B-field...



Direct transfer reactions possible in this region with 5 – 10 MeV/u beams from HIE-ISOLDE

The SpecMAT project:

To investigate exotic nuclei, specifically

- neutron-rich around $Z = 28$, e.g. $^{68, 70}\text{Ni}$; $^{78, 80}\text{Zn}$ (d, p); ($d, ^3\text{He}$)
- neutron-deficient around $Z = 82$, e.g. ^{184}Hg , ^{188}Pb , ^{196}Po (d, p); (p, d)
 - To place **active target** detector inside **solenoid magnet**,
 - & supplement with scintillator array in B-field...



Direct transfer reactions possible in this region with 5 – 10 MeV/u beams from HIE-ISOLDE

The SpecMAT project:

To investigate exotic nuclei, specifically

- neutron-rich around $Z = 28$, e.g. $^{68, 70}\text{Ni}$; $^{78, 80}\text{Zn}$ (d, p); ($d, ^3\text{He}$)
- neutron-deficient around $Z = 82$, e.g. ^{184}Hg , ^{188}Pb , ^{196}Po (d, p); (p, d)
 - To place active target detector inside solenoid magnet,
 - & supplement with **scintillator array** in B-field...



Direct transfer reactions possible in this region with 5 – 10 MeV/u beams from HIE-ISOLDE

LaBr₃(Ce) scintillators

Internal contamination problems due to Lanthanum

Density = **5.07 g/cm³**

Patented by Saint Gobain (though also sold through Canberra, Ortec..)

Price ~ **9 000 EUR** per scintillator for 1.5"*1.5" cylindrical model

Resolution = **2.7-3.3%** @ 662 keV

CeBr₃ scintillators


No internal background (apart from ²²⁷Ac contamination)

Density = **5.19 g/cm³**

Produced by Scionix Holland, RMD, Kinheng Crystal..

Price ~ **4 500 EUR** per scintillator for 1.5"*1.5" cylindrical model

Resolution = **3.8-4.0%** @ 662 keV



F.G.A. Quarati *et al.*, NIM A 729 (2013) 596-604



Cylindrical prototypes tested at Leuven with PMTs:

CeBr₃

- 1.5"×1.5"
- 2"×2"


LaBr₃

- 1.5"×1.5"



- Problem with scintillators + conventional PMTs:
- Detectors need to work in high B-field...

- Problem with scintillators + conventional PMTs:
- Detectors need to work in high B-field...
- Solution: replace PMTs with **Silicon Photomultipliers (SiPMs)**



Array-SMT User Manual – SensL, www.sensl.com/downloads/ds/UM-ArraySMT.pdf

LaBr₃(Ce) – tested at University of York

Industrial applications lab

Coupled 51x51x51 mm LaBr₃ to 8x8 array of 6 mm SiPMs

Optimized Voltage (29 to 30 V) and shaping time (0.5 μ s)



European Research Council
Established by the European Commission







SensL PCB





LaBr₃ crystal with
SiPMs




Crystal on PCB






2×2“ cubic crystal data obtained with SiPM and LaBr₃ (top),
and PM tube + LaBr₃ CeBr, NaI (bottom, F. Quarati *et al.* NIM A 729 (2013) 596-604)




- Good linearity of energy to signal





- Good linearity of energy to signal
- Resolution from fit, FWHM = 4.0% at 662 keV. After correction for non-linearity, FWHM = 4.2%

LaBr₃ with SiPMs vs PMTs




Tests to couple scintillators to GET electronics:

- Used Hamamatsu PMT, and
- Scintipack PMT base + HV source from ORTEC.
- Coupled to reduced CoBo setup at GANIL




Electronics



Preliminary data

Spectrum of ^{137}Cs with LaBr_3




Resolution of FWHM = 3.1% at 662 keV, out of nominal 2.6%

Preliminary data

Best resolution values with GET vs nominal values


Material	Dimensions	GET	Nominal
LaBr ₃	1.5"x1.5"	3.1%	2.6%
CeBr ₃	1.5"x1.5"	5.7%	4.0%
CeBr ₃	2"x2"	5.5%	4.0%

LaBr₃ + GET vs LaBr₃



Last test data to be obtained:

- Need crystals with **SiPMs + GET** electronics
- Workshop at Leuven designed adaptor PCB board for inner 6*6 pixels of 8*8 pixel SiPM board




Connections for read-out of a single pixel.
(Array-SMT User Manual – SensL,
www.sensl.com/downloads/ds/UM-ArraySMT.pdf)

- To test this setup with 1.5"*1.5" cubic LaBr₃ & CeBr₃
- Setup then comes to GANIL for tests with GET

Last test data to be obtained:

- Need crystals with **SiPMs + GET** electronics
- Workshop at Leuven designed adaptor PCB board for inner 6*6 pixels of 8*8 pixel SiPM board




Connections for read-out of a single pixel.
(Array-SMT User Manual – SensL,
www.sensl.com/downloads/ds/UM-ArraySMT.pdf)

- To test this setup with 1.5"*1.5" cubic LaBr₃ & CeBr₃
- Setup then comes to GANIL for tests with GET

Last test data to be obtained:

- Need crystals with **SiPMs + GET** electronics
- Workshop at Leuven designed adaptor PCB board for inner 6*6 pixels of 8*8 pixel SiPM board



Connections for read-out of a single pixel.
(Array-SMT User Manual – SensL,
www.sensl.com/downloads/ds/UM-ArraySMT.pdf)

- To test this setup with 1.5"*1.5" cubic LaBr₃ & CeBr₃
- Setup then comes to GANIL for tests with GET (circa Dec 2015)

SpecMAT – what is already measured:

- PMTs + Osprey module (LaBr_3 + CeBr_3)
- PMTs + GET (LaBr_3 + CeBr_3)
- SiPMs + analogue electronics (LaBr_3)

Not yet measured:

- SiPMs + GET
- Scintillators in B-field

Material	Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
LaBr3	1.5"×1.5"	cylinder		3.1%	2.8%	2.6%
LaBr3	2.0"×2.0"	cube	4.2%			
CeBr3	1.5"×1.5"	cylinder		5.7%	4.2%	4.0%
CeBr3	2.0"×2.0"	cylinder		5.5%	4.5%	4.0%

Measurements
Estimates

Measurements
Estimates

LaBr3

Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
1.5"×1.5"	cylinder		3,1%	2,8%	2,6%
2.0"×2.0"×2.0"	cube	4,2%			

Projection: SiPMs with GET	FWHM (at 662 keV)
4,7%	31
5,0%	33

CeBr3

Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
1.5"×1.5"	cylinder		5,7%	4,2%	4,0%
2.0"×2.0"	cylinder		5,5%	4,5%	4,0%
2.0"×2.0"×2.0"	cube - proj	5,2%			

Projection: SiPMs with GET	FWHM (at 662 keV)
7,1%	47
7,4%	49

LaBr3

Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
1.5"×1.5"	cylinder		3,1%	2,8%	2,6%
2.0"×2.0"×2.0"	cube	4,2%			

CeBr3

Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
1.5"×1.5"	cylinder		5,7%	4,2%	4,0%
2.0"×2.0"	cylinder		5,5%	4,5%	4,0%
2.0"×2.0"×2.0"	cube - proj	5,2%			

Projection: SiPMs with GET	FWHM (at 662 keV)
4,7 %	31
5,0 %	33
& in magnetic field:	
5,6 %	37
6,0 %	40

Projection: SiPMs with GET	FWHM (at 662 keV)
7,1 %	47
7,4 %	49
& in magnetic field:	
8,5 %	56
8,9 %	59

Measurements

Estimates

Measurements

Estimates

LaBr3

Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
1.5"×1.5"	cylinder		3,1%	2,8%	2,6%
2.0"×2.0"×2.0"	cube	4,2%			

CeBr3

Dimensions	Shape	SiPMs	GET	Osprey	Manufacturers
1.5"×1.5"	cylinder		5,7%	4,2%	4,0%
2.0"×2.0"	cylinder		5,5%	4,5%	4,0%
2.0"×2.0"×2.0"	cube - proj	5,2%			

Projection: SiPMs with GET	FWHM (at 662 keV)
4,7 %	31
5,0 %	33
& in magnetic field:	
5,6 %	37
6,0 %	40

Projection: SiPMs with GET	FWHM (at 662 keV)
7,1 %	47
7,4 %	49
& in magnetic field:	
8,5 %	56
8,9 %	59

Measurements

Estimates

Measurements

Estimates


- Two options: - Crystal + SiPM with mu-metal shield
- Crystal via light guide to PMT (outside solenoid)
- Doppler broadening also to be factored in

Acknowledgement:

The research leading to these results has received funding from the European Commission's Seventh Framework Programme (FP7/2007-2013) under the grant agreement SpecMAT (project no. 617156).

Efficiency measurements with cylindrical scintillators + PMTs



Efficiency vs E



Preliminary data

Best resolution values with GET vs nominal values

Material	Dimensions	GET	Nominal
LaBr ₃	1.5"x1.5"	3.1%	2.6%
CeBr ₃	1.5"x1.5"	5.7%	4.0%
CeBr ₃	2"x2"	5.5%	4.0%



Efficiency measurements with cylindrical scintillators + PMTs

